

Bruno Scooter-Lift II

Bruno Model VSL-900

Lift and store your fully assembled scooter or wheelchair through either side sliding door of minivans.

Only Bruno Independent Living Aids, Inc.[®] has over 18 different scooter, wheelchair or power chair lifts for nearly any car, van, truck or sport utility vehicle.

Plus, all Bruno lifts are eligible for reimbursement under the terms of Mobility Programs offered by GM, Saturn, Ford, and Chrysler.



- ◆ The Scooter-Lift II can raise and store a fully assembled scooter or wheelchair through the driver's or passenger's side sliding door of minivans, and passenger side doors of full size vans.
- ◆ The Scooter-Lift II requires very little arm strength, balance or effort to use. It is so compact that even in a minivan, there is still enough room for three passengers to sit comfortably in the rear seat. With scooter removed, or when travelling without the scooter, entry and exit is easy since the door opening remains unobstructed.
- ◆ Bruno's Scooter-Lift II also features an optional installation assembly that eliminates drilling any holes in the floor of the vehicle. This means easier removal when it comes time to transfer the lift in and out of a minivan. Contact **Bruno's Technical Applications Department** about specific applications for this system.
- ◆ The Scooter-Lift II can be used for either three- or four-wheel scooters and manual folding or power wheelchairs weighing up to 200 lbs. (91k). The Scooter-Lift II's hand-held controller makes this lift simple and easy to use. The Scooter-Lift II uses electrical power from the vehicle's battery and all lift operations are automatic. You only need to help guide the scooter or wheelchair into place.
- ◆ Expert installation and complete after-sale service is provided by Bruno's nationwide and international network of dealers.



Bruno Independent Living Aids, Inc., 1780 Executive Drive
P.O. Box 84, Oconomowoc, WI 53066

1-262-567-4990 ◆ FAX: 1-262-953-5501 ◆ www.bruno.com

Using the *Bruno Scooter-Lift II* is easy!



1 Push the **OUT** button to bring the lift above the scooter at the side door, then hook the C-arm to the T-bar docking device.



2 To lift the scooter, push the **UP** button until the wheels clear the van's floor.



3 Once the scooter is at the proper height, push the **IN** button.



4 When the scooter is in the van, push the **DOWN** button to set it on the floor.



5 The Scooter-Lift II also powers a fully assembled wheelchair into the van effortlessly.



6 Bruno manufactures a version of the Scooter-Lift II for virtually all minivans with sliding doors, and also for the passenger's side of full-size vans.

Represented By:



Bruno Independent Living Aids, Inc., 1780 Executive Drive,
P.O. Box 84, Oconomowoc, WI 53066

1-262-567-4990 ♦ FAX: 1-262-953-5501
www.bruno.com